



EGIS Non-Vascular Stent

Airway Stent

Stents for trachea and bronchus

The Airway Stent

EGIS Airway Stent is used for maintaining the patency of the trachea and the Bronchus. The two points, proximal and distal GOLD MARKERS, facilitate the procedure and visibility of the stent under fluoroscopy.

Airway Stent is used for the treatment of airway strictures and fistulas, caused by malignant obstruction or congenital deformities. As a result, patients will be able to breathe normally. The stent is designed to have excellent radial force against rigid trachea and bronchial strictures.

Specification

Ordering No. First two** : Diameter(mm) / Last two** : Length(cm)

Ordering Information		Size	
Type	Ordering No.	Diameter	Length
BRONCHUS	Straight Bare	10, 12, 14	3, 4, 5, 6
	Straight Covered		
	Single Flare Covered		
	Double Flare Covered		
TRACHEA	Straight Covered	16, 18, 20	4, 5, 6, 7, 8, 9
	Single Flare Covered		

