



Non–Vascular Stent

Colorectal/Pyloric Stent

The Colorectal Stent

EGIS Colorectal Stent is used to open strictures in the large intestine. The Stent is designed with outstanding conformability to fit into the large intestine smoothly. Colorectal Stent naturally fits into the lower digestive tract. It allows patients to defecate normally. The proximal, middle and distal GOLD markers increase visibility under fluoroscopy.

Delivery System 10Fr 120/180cm

Specification

Ordering No. First two** : Diameter(mm) / Last two** : Length(cm)

Type	Ordering No.		Size	
	OTW	TTS	Diameter	Length
Single Bare	CSB****08	CSB****22	18, 20, 22, 24 26, 28 30	6, 8, 10, 12, 14 (*Only Bare)
	CSB****15			
Double Bare	CDB****08	CDB****22		
	CDB****15			
Single Covered	CSC****08	CSC****22		
	CSC****15			
Double Covered	CDC****08	CDC****22		
	CDC****15			

The Pyloric Stent

EGIS Pyloric Stent is designed to fit into the gastric outlet and duodenum naturally. The proximal, middle and distal GOLD markers increase visibility under fluoroscopy. The optimal radial force of the EGIS Pyloric stent will allow patients to achieve normal digestion. The stent is designed for maintaining excellent conformability in order to adjust to the shape of the gastric outlet and duodenum smoothly.

Delivery System 10Fr, 80/150/220cm, 12Fr, 80/150/220cm

Specification

Ordering No. First two** : Diameter(mm) / Last two** : Length(cm)

Type	Ordering No.		Size	
	OTW	TTS	Diameter	Length
Single Bare	PSB****12	PSB****18	18, 20, 22, 24	6, 8, 10, 12 14 (*Only Bare)
Double Bare	PDB****12	PDB****18		
Single Covered	PSC****12	PSC****18		
Double Covered	PDC****12	PDC****18		